

Protocols - ^{13}C Breath Tests - STOMACH




Table of Contents

<i>Introduction</i>	3
• ¹³ C Breath Test protocols	3
• ¹³ C Breath Testing: principle and requirements.....	3
• Preparation of the patient.....	3
• Administration of ¹³ C labeled substrate	3
• Collection of breath samples	3
• Measurement of ¹³ C enrichment.....	3
• Calculation of the end result.....	3
• Applications	4
• Literature	4
• Note.....	4
• Authors	4
 <i>Test Protocols Stomach</i>	 5
• Helicobacter pylori Infection	5
¹³ C-Urea Breath Test	
• Gastric Emptying of Solids	6
¹³ C-Octanoic Acid Breath Test (Leuven Model)	
• Gastric Emptying of Liquids.....	8
¹³ C-Acetate Breath Test	
 <i>Literature</i>	 9
Stomach	9
• Helicobacter pylori Infection	9
¹³ C-Urea Breath Test	
• Gastric Emptying of Solids	13
¹³ C-Octanoic Acid Breath Test (Leuven Model)	
• Gastric Emptying of Liquids.....	15
¹³ C-Acetate Breath Test	

Introduction

■ **¹³C Breath Test protocols**

This folder contains a set of protocols describing the principles and general test procedures for today's most relevant ¹³C Breath Tests to study specific functions of the liver, pancreas, stomach and intestine. The list will be updated regularly adding additional tests or additional information on the already described tests.

The information is meant as a start to enter the field of stable isotope ¹³C Breath testing initiated by the interest in a specific test. To actually be able to introduce a test in your hospital you must familiarize yourself with basic knowledge of breath testing with ¹³C substrates and the existing knowledge on the particular application of interest. There is no such thing as a standard protocol for all tests.

■ **¹³C Breath Testing: principle and requirements**

A ¹³C Breath Test consists of the administration to a patient of a ¹³C labeled substrate that is metabolized by a specific enzyme system resulting in ¹³CO₂ as the end product. To monitor the enzyme response ¹³C enrichment is measured in breath CO₂.

The total procedure of ¹³C Breath testing includes the definition of the preparation of the patient before the test, administration of the ¹³C labeled substrate, collection of breath samples, measurement of ¹³C enrichment in breath CO₂ and the calculation of an end result.

■ **Preparation of the patient**

In general, tests will be performed in the fasting state and the patient should be at a low and stable level of natural ¹³C abundance. Therefore, the patient must be instructed to avoid eating ¹³C enriched foods such as corn products, cane sugar, pineapple and tequila the last days before the test and to come to the clinic fasted. In certain cases (¹³C Lactose-Ureide breath test) the patient must be pretreated with unlabeled substrate to stimulate the involved enzyme system.

■ **Administration of ¹³C labeled substrate**

The test substrate may be administered as a simple solution in water with or without a standardized test meal. Sometimes it needs to be incorporated into a specific ingredient of the meal. The test meal and the dose of substrate may be different for adults and children.

■ **Collection of breath samples**

Every protocol has its own time schedule of breath collections. The number of samples may be as small as 2 or more than 20. To define the ¹³C enrichment in breath CO₂, it is also necessary to obtain at least two breath samples before the ingestion of the ¹³C substrate to determine the natural background of ¹³C abundance. The methodology of collecting breath samples is dependent on the technology to determine the ¹³C enrichment. The protocols are based on Continuous Flow Isotope Ratio Mass Spectrometry as the analytical technique. In this case breath samples are simply blown through a straw into special 10 ml gas collection tubes that directly fit into the sample tray of the instruments. In case of Infrared technology, special bags provided by the instrument manufacturer must be used.

■ **Measurement of ¹³C enrichment**

To determine the ¹³C abundance in breath CO₂ you need the availability of Isotope Ratio Mass Spectrometry (IRMS) or specialized Infrared instrumentation. The protocols are based on Isotope Ratio Mass Spectrometry. For a number of tests (Aminopyrine, Methacetin, Urea) Infrared Spectroscopy has proven to be a valid alternative analytical technique. For other tests Infrared technology has not yet been validated so far. In principle, the test substrate is not a determinant of the validity of the analytical technique. It is the level of ¹³C enrichment that determines the analytical requirement. Validation of Infrared analysis for other application is recommended, as it is recommended to validate any breath test in your own clinical laboratory. You may have instrumentation available or contact a service center for the analyses.

■ **Calculation of the end result**

For some tests the only calculation needed is the subtraction of the natural background value from the measured value at a defined time. In other cases it is necessary to calculate the amount of ¹³C that is recovered in breath during the experimental period. In a third type of application the time course of the enrichment appearance is of importance requiring calculation of the rate of appearance.

■ Applications

In the present update the following tests have been described:

Table of contents of applications

Function		¹³ C Substrate	Page
Stomach (Registered)			
1.	Helico pylori Infection	¹³ C-Urea	5
Stomach			
2.	Gastric Emptying of Solids (Leuven Model)	¹³ C-Octanoic Acid	6
3.	Gastric Emptying of Liquids	¹³ C-Sodium Acetate	8
Literature			
			9

■ Literature

Included is a list of literature references that will introduce you to the most important articles describing aspects of the different tests described in the protocols.

■ Note

Great care has been taken over the composition of the text, figures and tables. The possibility of errors however, cannot be excluded completely. Therefore Campro Scientific GmbH and the authors cannot accept any legal or other liability with respect to incorrect details and their consequences. The authors would be grateful to receive suggestions for improvements and information about errors. This information may not be photocopied, duplicated or translated into another language neither partially nor completely, in any form without the written agreement of Campro Scientific GmbH.

■ Authors

1. Dr. F. Stellaard
University Hospital Groningen, the Netherlands
Dept. Pathology and Laboratory Medicine
And Centre for Liver, Intestinal and Metabolic Disease
2. Dr. Ahmad Rajabi
Campro Scientific GmbH
Berlin, Germany

© Campro Scientific GmbH

European Headquarters

P.O. Box 37 03 31
D-14133 Berlin
Germany
Tel. : +49(0)30.629.01.89.0
Fax : +49(0)30.629.01.89.89
E-mail : info@campro.eu
Web : www.campro.eu

Dutch Sales Office

P.O. Box 316
NL-3900 AH Veenendaal
The Netherlands
Tel. : +31(0)318.529.437
Fax : +31(0)318.542.181
E-mail : info@campro.eu
Web : www.campro.eu

Please contact us for technical and price information.

Test Protocols Stomach

1. Helicobacter pylori Infection

¹³C-Urea Breath Test

■ Principle

¹³C-Urea contains one carbon atom labeled with the non-radioactive isotope ¹³C. After oral administration, ¹³C-Urea passes through the stomach. The presence of Helicobacter pylori leads to the hydrolysis of ¹³C-Urea to ammonia and ¹³CO₂ due to the action of urease enzyme. ¹³CO₂ is absorbed, transported by the blood and excreted by the lungs. The appearance of ¹³C in breath CO₂ reflects the presence of Helicobacter pylori.

■ Applicability of ¹³C-Urea Breath Test

¹³C-Urea Breath Test has so far been applied to adults and children.

■ Applications

¹³C-Urea Breath Test is used to detect Helicobacter pylori in patients with general gastric complaints. An increase of ¹³C abundance outside the control range is indicative for a positive test result. The effect of eradication therapy can be monitored within one month after the completion of therapy.

■ Protocol

Adults: The ¹³C-Urea Breath Test is performed after an overnight fast. A dose of 75 mg or 100 mg is administered orally after dissolution in about 100 ml water. The intake of the solution is preceded by citric acid solution or orange juice. Breath samples are collected before (2x) and 30 minutes (2x) after ingestion of the ¹³C-Urea. ¹³C enrichment in breath CO₂ is determined by Isotope Ratio Mass Spectrometry (IRMS) or non-dispersive Infrared Spectrometry. The absolute increase of ¹³C abundance after 30 minutes is compared to the baseline value and then used as the diagnostic parameter.

Children: A dose of 50 mg is sufficient for children. The same time schedule for breath collections can be used.

■ Interpretation of test results

It is advised to obtain your own internal control values. Generally a cut-off value of 4 or 5 ‰ is described for the increase after 30 minutes over the baseline value before administration of ¹³C-Urea.

■ Precautions

The results must be interpreted with caution when the patient has undergone gastrectomy, has been treated with antibiotics, proton pump inhibitors and antisecretory drugs, or in other cases in which acid secretion is affected such as atrophic gastritis.

■ Summary

	Dose	Samples	
Adults	75 mg ¹³ C-Urea or 100 mg ¹³ C-Urea	2	Before administration
		2	30 minutes (0.5 h) after administration
Children	50 mg ¹³ C-Urea	2	Before administration
		2	30 minutes (0.5 h) after administration

2. Gastric Emptying of Solids

¹³C-Octanoic Acid Breath Test (Leuven Model)

■ Principle

(1-¹³C)-Octanoic Acid contains a carboxyl-carbon labeled with the non-radioactive isotope ¹³C. After oral administration (1-¹³C)-Octanoic Acid passes through the stomach and is completely absorbed by the small intestine. (1-¹³C)-Octanoic Acid is oxidized to a large extent. The kinetics of appearance of ¹³C in breath CO₂ reflects the rate of gastric emptying of the solid phase of a meal.

■ Applicability of ¹³C-Octanoic Acid Breath Test

¹³C-Octanoic Acid Breath Test has so far been applied to adults, children and newborn infants.

■ Applications

¹³C-Octanoic Acid Breath test is used to detect impaired gastric emptying of solids by comparison with a control range. The effect of drug treatment can be monitored.

■ Protocol

Adults: The ¹³C-Octanoic Acid Breath Test is performed after an overnight fast. A dose of 100 mg (1-¹³C)-Octanoic Acid is administered orally in a solid test meal. The test meal is standardized and consists of one scrambled egg with two slices of white bread and 5 g of margarine, together with 150 ml water (swallowed immediately after ingestion of the meal). The total caloric content is 250 kcal. The egg yolk is doped with 100 mg (1-¹³C)-Octanoic Acid and fried separately from the egg white. The meal is consumed within 10 minutes. Breath samples are collected before (2x), every 5 minutes during the first 30 minutes (0.5 h) and every 15 minutes for the next 210 minutes (3.5 h) after the ingestion of the (1-¹³C)-Octanoic Acid. ¹³C enrichment in breath CO₂ is determined by Isotope Ratio Mass Spectrometry (IRMS). The equation of the breath test results is obtained by 2 non-linear regression curves fitting the % dose ¹³C recovered in breath per minute and the cumulative % dose recovered in breath. From this equation the half emptying time and the lag phase time are calculated as well as the gastric emptying coefficient (GEC).

Children: A pancake is made of 5 g sugar, 12.5 g flour, 10 g full cream milk and one egg. Half of the egg yolk is baked separately after solubilising 50 mg (1-¹³C)-Octanoic Acid by mixing. The obtained paste is baked around the labeled half of the egg yolk. Total energy content is 150 kcal. Meal consumption, breath sampling and calculations are the same as described above for adults.

Newborn infants: A volume of water, necessary for the baby's usual intake, is heated. 50 mg (1-¹³C)-Octanoic Acid and 1 g PolyEthylene Glycol 3350 needed for solubilization are mixed and cooked. The adequate amount of milk powder is added, and the bottle kept in the refrigerator overnight. The infants receive the test meal at the time of their usual first morning feeding, which allows for 3 hours fasting. The test meal is re-cooked and cooled prior to administration. Breath samples are collected using a nasal prong carefully inserted into the nasopharynx. 10 ml exhaled air is slowly aspirated in a syringe during the expiration phase. Breath samples are collected before (2x), every 5 minutes during the first 30 minutes (0.5 h) and every 15 minutes during the next 210 minutes (3.5 h) after drinking the meal.

■ Interpretation of test results

It is advised to obtain your own internal control values. The control values will depend on the population (age group) and the test meal used.

■ Precautions

No contra-indications for the ¹³C-Octanoic Acid Breath Test have been described so far.

■ Summary

	Dose	Samples	
Adults	100 mg (1- ¹³ C)-Octanoic Acid	2	Before administration
		6	Every 5 minutes for the first 30 minutes after administration (0.5.h)
		14	Every 15 minutes for the next 210 minutes after administration (3.5 h)
Children / Newborn	50 mg (1- ¹³ C)-Octanoic Acid	2	Before administration
		6	Every 5 minutes for the first 30 minutes after administration (0.5.h)
		14	Every 15 minutes for the next 210 minutes after administration (3.5 h)

3. Gastric Emptying of Liquids

¹³C-Acetate Breath Test

■ Principle

(1-¹³C)-Sodium Acetate contains a carboxyl-carbon labeled with the non-radioactive isotope ¹³C. After oral administration (1-¹³C)-Sodium Acetate passes through the stomach and is completely absorbed by the small intestine. (1-¹³C)-Sodium Acetate is oxidized to a large extent. The kinetics of appearance of ¹³C in breath CO₂ reflects the rate of gastric emptying of the liquid phase of a meal.

■ Applicability of ¹³C-Acetate Breath Test

¹³C-Acetate Breath Test has so far been applied to adults and children.

■ Applications

¹³C-Acetate Breath Test is used to detect impaired gastric emptying of liquids by comparison with a control range. The effect of drug treatment can be monitored.

■ Protocol

The ¹³C-Acetate Breath Test is performed after an overnight fast. A dose of 150 mg (1-¹³C)-Sodium Acetate is administered orally in a liquid test meal. So far the test meal has not yet been standardized and should be adapted to the individual question (neonatal feeding, exercise experiments, food matrix studies). Breath samples are collected before (2x) and every 5 minutes for the next 120 minutes (2 h) and every 10 minutes for the following 120 minutes (2 h) after the ingestion of the (1-¹³C)-Sodium Acetate. ¹³C enrichment in breath CO₂ is determined by Isotope Ratio Mass Spectrometry (IRMS). The equation of the breath test results is obtained by non-linear regression curve, fitting of the % dose ¹³C recovered in breath per minute. From this equation the half emptying time and the lag phase are calculated. Due to a high correlation between the half emptying time and the time of the peak maximum of ¹³C appearance in breath, the time to peak maximum appears a simple diagnostic parameter for the interpretation of the rate of gastric emptying.

■ Interpretation of test results

It is advised to obtain your own internal control values. The control values will depend on the population (age group) and the test meal that has been used. For accurate kinetic data the half emptying rate should be calculated. For screening purposes, the time to peak maximum is a simple but accurate marker.

■ Precautions

No contra-indications for the ¹³C-Acetate Breath Test test have been described so far.

■ Summary

	Dose	Samples	
	150 mg (1- ¹³ C)- Sodium Acetate	2	Before administration
		24	Every 5 minutes, the first 120 minutes after administration (2 h)
		12	Every 10 minutes, the second 120 minutes after administration (2 h)

Literature

Stomach

1. Helicobacter pylori Infection ¹³C-Urea Breath Test

■ Recommended literature

Recent literature

1. Dulbecco P, Gambaro C, Bilardi C, Zentilin P, Mele MR, Mansi C, Biagini R, Tessieri L, Iiritano E, Usai P, Vigneri S, Savarino V. Impact of long-term ranitidine and pantoprazole on accuracy of [¹³C]urea breath. test. Dig Dis Sci. 2003 Feb;48(2):315-21.
2. Gatta L, Vakil N, Ricci C, Osborn JF, Tampieri A, Perna F, Miglioli M, Vaira D. A rapid, low-dose, ¹³C-urea tablet for the detection of Helicobacter pylori infection before and after treatment. Aliment Pharmacol Ther. 2003 Mar;17(6):793-8.
3. Schmidt G. Carbon isotope analysis in urea at high ¹³C-abundances using the ¹³/¹²CO₂-breath test device FANci2. Isotopes Environ Health Stud. 2002 Sep;38(3):185-8.
4. Opekun AR, Abdalla N, Sutton FM, Hammoud F, Kuo GM, Torres E, Steinbauer J, Graham DY. Urea breath testing and analysis in the primary care office. J Fam Pract. 2002 Dec;51(12):1030-2.
5. Wong WM, Lam SK, Lai KC, Chu KM, Xia HH, Wong KW, Cheung KL, Lin SK, Wong BC. A rapid-release 50-mg tablet-based ¹³C-urea breath test for the diagnosis of Helicobacter pylori infection. Aliment Pharmacol Ther. 2003 Jan;17(2):253-7.
6. Saltik IN, Demir H, Kocak N, Ozen H, Gurakan F, Yuce A. Diagnostic accuracy of ¹³C-urea breath test for Turkish children with Helicobacter pylori infection. Am J Gastroenterol. 2003 Jan;98(1):222-3.
7. Chua TS, Fock KM, Teo EK, Ng TM. Validation of ¹³C-urea breath test for the diagnosis of Helicobacter pylori infection in the Singapore population. Singapore Med J. 2002 Aug;43(8):408-11.
8. Canete A, Abunaji Y, Alvarez-Calatayud G, DeVicente M, Gonzalez-Holguera JA, Leralta M, Pajares JM, Gisbert JP. Breath Test Using A Single 50-mg Dose of ¹³C-Urea to Detect Helicobacter pylori Infection in Children. J Pediatr Gastroenterol Nutr. 2003 Jan;36(1):105-11.
9. Perri F, Manes G, Neri M, Vaira D, Nardone G. Helicobacter pylori antigen stool test and ¹³C-urea breath test in patients after eradication treatments. Am J Gastroenterol. 2002 Nov;97(11):2756-62.
10. Kurpad AV, Ajami A, Young VR. ¹³C breath tests in infections and beyond. Food Nutr Bull. 2002 Sep;23(3 Suppl):21-9.
11. Isomoto H, Inoue K, Shikuwa S, Furusu H, Nishiyama T, Omagari K, Mizuta Y, Murase K, Murata I, Enjoji A, Kanematsu T, Kohno S. Five minute endoscopic urea breath test with 25 mg of (¹³C)-urea in the management of Helicobacter pylori infection. Eur J Gastroenterol Hepatol. 2002 Oct;14(10):1093-100.
12. Gee I, Playford RJ, Turner D, Sheldon N, Wicks AC. Cost analysis of breath test versus endoscopy for dyspepsia. Digestion. 2002;65(4):207-12.

13. Perri F, Ricciardi R, Merla A, Piepoli A, Gasperi V, Quitadamo M, Andriulli A. Appropriateness of urea breath test: a prospective observational study based on Maastricht 2000 guidelines. *Aliment Pharmacol Ther.* 2002 Aug;16(8):1443-7.
14. Leung WK, Hung LC, Kwok CK, Leong RW, Ng DK, Sung JJ. Follow up of serial urea breath test results in patients after consumption of antibiotics for non-gastric infections. *World J Gastroenterol.* 2002 Aug;8(4):703-6.
15. Kato S, Ozawa K, Konno M, Tajiri H, Yoshimura N, Shimizu T, Fujisawa T, Abukawa D, Minoura T, Iinuma K. Diagnostic accuracy of the ¹³C-urea breath test for childhood *Helicobacter pylori* infection: a multicenter Japanese study. *Am J Gastroenterol.* 2002 Jul;97(7):1668-73.
16. Kubota K, Shimoyama S, Shimizu N, Noguchi C, Mafune K, Kaminishi M, Tange T. Studies of ¹³C-urea breath test for diagnosis of *Helicobacter pylori* infection in patients after partial gastrectomy. *Digestion.* 2002;65(2):82-6.
17. Herold R, Becker M. ¹³C-urea breath test threshold calculation and evaluation for the detection of *Helicobacter pylori* infection in children. *BMC Gastroenterol.* 2002 May 6;2(1):12.
18. Franceschi F, Armuzzi A, Cremonini F, Carloni E, Zocco MA, Di Caro S, Padalino C, Genta RM, Pola P, Gasbarrini G, Gasbarrini A. Delta¹³CO₂ excretion and expression of dyspeptic symptoms in patients evaluated for *Helicobacter pylori* infection by [¹³C] urea breath test. *Dig Dis Sci.* 2002 Apr;47(4):804-8.
19. Hegedus O, Ryden J, Rehnberg AS, Nilsson S, Hellstrom PM. Validated accuracy of a novel urea breath test for rapid *Helicobacter pylori* detection and in-office analysis. *Eur J Gastroenterol Hepatol.* 2002 May;14(5):513-20.
20. Bode G, Hoffmeister A, Koenig W, Brenner H, Rothenbacher D. Characteristics of differences in *Helicobacter pylori* serology and ¹³C-urea breath-testing in an asymptomatic sample of blood donors. *Scand J Clin Lab Invest.* 2001;61(8):603-8.
21. Yoshimura N, Tajiri H, Sawada A, Kozaiwa K, Ida S, Fujisawa T, Konno M, Kato S. A ¹³C-urea breath test in children with *Helicobacter pylori* infection: assessment of eradication therapy and follow-up after treatment. *J Gastroenterol.* 2001 Sep;36(9):606-11.
22. Chen X, Haruma K, Kamada T, Hartori N, Yoshihara M, Kitadai Y, Tanaka S, Sumii K, Chayama K. A low ¹³C-urea breath test value is associated with increased risk of gastric cancer. *J Gastroenterol.* 2001 Sep;36(9):601-5.
23. Yeh JL, Peng YC, Tung CF, Chen GH, Chow WK, Chang CS, Yeh HZ, Poon SK. Role of *Helicobacter pylori* in cirrhotic patients with dyspepsia: a ¹³C-urea breath test study. *Adv Ther.* 2001 May-Jun;18(3):140-50.
24. Wong WM, Wong BC, Li TM, Wong KW, Cheung KL, Fung FM, Xia HH, Lam SK. Twenty-minute 50 mg ¹³C-urea breath test without test meal for the diagnosis of *Helicobacter pylori* infection in Chinese. *Aliment Pharmacol Ther.* 2001 Sep;15(9):1499-504.
25. Leodolter A, Wolle K, Malfertheiner P. Current standards in the diagnosis of *Helicobacter pylori* infection. *Dig Dis.* 2001;19(2):116-22. Review.

26. Gisbert JP, Vazquez MA, Jimenez I, Cruzado AI, Carpio D, Del Castillo E, Martin MJ, Morales A, Pajares R, Rodriguez A, Pajares JM. ¹³C-urea breath test for the diagnosis of *Helicobacter pylori* infection before treatment: is citric acid necessary? *Dig Liver Dis*. 2000 Jan-Feb;32(1):20-4.
27. Malaty HM, Logan ND, Graham DY, Ramchatesingh JE, Reddy SG. *Helicobacter pylori* infection in asymptomatic children: comparison of diagnostic tests. *Helicobacter*. 2000 Sep;5(3):155-9.
28. Sheu BS, Lee SC, Yang HB, Kuo AW, Wang YL, Shiesh SC, Wu JJ, Lin XZ. Selection of lower cutoff point of [¹³C]urea breath test is helpful to monitor *H. pylori* eradication after proton pump inhibitor-based triple therapy. *Dig Dis Sci*. 2000 Jul;45(7):1330-6.
29. Suto G, Vincze A, Pakodi F, Hunyady B, Karadi O, Garamszegi M, Laszlo T, Mozsik G. ¹³C-Urea breath test is superior in sensitivity to detect *Helicobacter pylori* infection than either antral histology or rapid urease test. *J Physiol Paris*. 2000 Mar- Apr;94(2):153-6.
30. Peng NJ, Hsu PI, Lee SC, Tseng HH, Huang WK, Tsay DG, Ger LP, Lo GH, Lin CK, Tsai CC, Lai KH. A 15-minute [¹³C]-urea breath test for the diagnosis of *Helicobacter pylori* infection in patients with non-ulcer dyspepsia. *J Gastroenterol Hepatol*. 2000 Mar;15(3): 284-9.
31. Bazzoli F, Cecchini L, Corvaglia L, Dall'Antonia M, De Giacomo C, Fossi S, Casali LG, Gullini S, Lazzari R, Leggeri G, Lerro P, Valdambrini V, Mandrioli G, Marani M, Martelli P, Miano A, Nicolini G, Oderda G, Pazzi P, Pozzato P, Ricciardiello L, Roda E, Simoni P, Sottili S, Zagari RM. Validation of the ¹³C-urea breath test for the diagnosis of *Helicobacter pylori* infection in children: a multicenter study. *Am J Gastroenterol*. 2000 Mar;95(3):646-50.
32. Riepl RL, Folwaczny C, Otto B, Klauser A, Blendinger C, Wiebecke B, Konig A, Lehnert P, Heldwein W. Accuracy of ¹³C-urea breath test in clinical use for diagnosis of *Helicobacter pylori* infection. *Z Gastroenterol*. 2000 Jan;38(1):13-9.
33. Kindermann A, Demmelmair H, Koletzko B, Krauss-Etschmann S, Wiebecke B, Koletzko S. Influence of age on ¹³C-urea breath test results in children. *J Pediatr Gastroenterol Nutr*. 2000 Jan;30(1):85-91.
34. Vincent P, Michaud L, Martin de Lasalle E, Benon B, Turck D, Gottrand F. ¹³C-urea breath test and gastric mucosal colonization by *Helicobacter pylori* in children: quantitative relation and usefulness for diagnosis of infection. *Helicobacter*. 1999 Dec;4(4):233-7.
35. Corvaglia L, Bontems P, Devaster JM, Heimann P, Glupczynski Y, Keppens E, Cadranet S. Accuracy of serology and ¹³C-urea breath test for detection of *Helicobacter pylori* in children. *Pediatr Infect Dis J*. 1999 Nov;18(11):976-9.
36. Connor SJ, Seow F, Ngu MC, Katelaris PH. The effect of dosing with omeprazole on the accuracy of the ¹³C-urea breath test in *Helicobacter pylori*-infected subjects. *Aliment Pharmacol Ther*. 1999 Oct;13(10):1287-93.
37. Connor SJ, Ngu MC, Katelaris PH. The impact of short-term ranitidine use on the precision of the ¹³C-urea breath test in subjects infected with *Helicobacter pylori*. *Eur J Gastroenterol Hepatol*. 1999 Oct;11(10):1135-8.

38. Colaiocco Ferrante L, Papponetti M, Marcuccitti J, Neri M, Festi D. ¹³C-urea breath test for helicobacter pylori infection: stability of samples over time. Scand J Gastroenterol. 1999 Sep;34(9):942-3.
39. Bravo LE, Realpe JL, Campo C, Mera R, Correa P. Effects of acid suppression and bismuth medications on the performance of diagnostic tests for Helicobacter pylori infection. Am J Gastroenterol. 1999 Sep;94(9):2380-3.
40. Zagari RM, Bazzoli F, Pozzato P, Fossi S, De Luca L, Nicolini G, Berretti D, Roda E. Review article: non-invasive methods for the diagnosis of Helicobacter pylori infection. Ital J Gastroenterol Hepatol. 1999 Jun-Jul;31(5):408-15.

2. Gastric Emptying of Solids

¹³C-Octanoic Acid Breath Test

■ Recommended literature

1. Bromer MQ, Kantor SB, Wagner DA, Knight LC, Maurer AH, Parkman HP. Simultaneous measurement of gastric emptying with a simple muffin meal using [¹³C]octanoate breath test and scintigraphy in normal subjects and patients with dyspeptic symptoms. *Dig Dis Sci.* 2002 Jul;47(7): 1657-63.
2. Gonlachanvit S, Chey WD, Goodman KJ, Parkman HP. Effect of meal size and test duration on gastric emptying and gastric myoelectrical activity as determined with simultaneous [¹³C]octanoate breath test and electrogastrography in normal subjects using a muffin meal. *Dig Dis Sci.* 2001 Dec;46(12):2643-50.
3. Ritz MA, Fraser R, Edwards N, Di Matteo AC, Chapman M, Butler R, Cmielewski P, Tournadre JP, Davidson G, Dent J. Delayed gastric emptying in ventilated critically ill patients: measurement by ¹³C-octanoic acid breath test. *Crit Care Med.* 2001 Sep;29(9):1744-9.
4. Mion F, Ecochard R, Guitton J, Ponchon T. (¹³CO₂) breath tests: comparison of isotope ratio mass spectrometry and non-dispersive infrared spectrometry results. *Gastroenterol Clin Biol.* 2001 Apr;25(4):375-9.
5. Kulik W, van Weissenbruch MM, Menelik N, Cranendonk A, Kneepkens CM, Lafeber HN. Improved use of the [¹³C]octanoic acid breath test as intra-individual parameter to study the effect of a prokinetic drug on gastric emptying in preterm infants with oral feeding intolerance. *J Chromatogr B Biomed Sci Appl.* 2001 Jan 5;750(1):147-53.
6. Chey WD, Shapiro B, Zawadski A, Goodman K. Gastric emptying characteristics of a novel (¹³C)-octanoate-labeled muffin meal. *J Clin Gastroenterol.* 2001 May-Jun;32(5): 394-9.
7. Kim DY, Myung SJ, Camilleri M. Novel testing of human gastric motor and sensory functions: rationale, methods, and potential applications in clinical practice. *Am J Gastroenterol.* 2000 Dec;95(12):3365-73. Review.
8. Cappello G, Malatesta MG, Ferri A, Ciccaglione AF, Toracchio S, Grossi L, Marzio L. Gastric emptying of a solid-liquid meal measured with ¹³C octanoic acid breath test and real-time ultrasonography: a comparative study. *Am J Gastroenterol.* 2000 Nov;95(11):3097-100.
9. Peracchi M, Gebbia C, Ogliari C, Fraquelli M, Viganò R, Baldassarri A, Bianchi PA, Conte D. Influence of caloric intake on gastric emptying of solids assessed by ¹³C-octanoic acid breath test. *Scand J Gastroenterol.* 2000 Aug;35(8):814-8.
10. Robertson MD, Mathers JC. Gastric emptying rate of solids is reduced in a group of ileostomy patients. *Dig Dis Sci.* 2000 Jul;45(7):1285-92.
11. Delbende B, Perri F, Couturier O, Leodolter A, Mauger P, Bridgi B, Bizais Y, des Varannes SB, Andriulli A, Galmiche JP. ¹³C-octanoic acid breath test for gastric emptying measurement. *Eur J Gastroenterol Hepatol.* 2000 Jan;12(1):85-91.
12. Van Den Driessche M, Peeters K, Marien P, Ghooys Y, Devlieger H, Veereman-Wauters G. Gastric emptying in formula-fed and breast-fed infants measured with the ¹³C-octanoic acid breath test. *J Pediatr Gastroenterol Nutr.* 1999 Jul;29(1):46-51.
13. Barnett C, Snel A, Omari T, Davidson G, Haslam R, Butler R. Reproducibility of the ¹³C-octanoic acid breath test for assessment of gastric emptying in healthy preterm infants. *J Pediatr Gastroenterol Nutr.* 1999 Jul;29(1):26-30.

14. Perri F, Clemente R, Festa V, Quitadamo M, Niro G, Andriulli A. ¹³C-octanoic acid breath test: a reliable tool for measuring gastric emptying. *Ital J Gastroenterol Hepatol.* 1998 Apr;30(2):211-7. Review.
15. Verhagen MA, Samsom M, Maes B, Geypens BJ, Ghoo YF, Smout AJ. Effects of a new motilide, ABT-229, on gastric emptying and postprandial antroduodenal motility in healthy volunteers. *Aliment Pharmacol Ther.* 1997 Dec;11(6):1077-86.
16. Choi MG, Camilleri M, Burton DD, Zinsmeister AR, Forstrom LA, Nair KS. Reproducibility and simplification of ¹³C-octanoic acid breath test for gastric emptying of solids. *Am J Gastroenterol.* 1998 Jan;93(1):92-8.
17. Choi MG, Camilleri M, Burton DD, Zinsmeister AR, Forstrom LA, Nair KS. [¹³C]octanoic acid breath test for gastric emptying of solids: accuracy, reproducibility, and comparison with scintigraphy. *Gastroenterology.* 1997 Apr;112(4):1155-62.
18. Veereman-Wauters G, Ghoo Y, van der Schoor S, Maes B, Hebbalkar N, Devlieger H, Eggermont E. The ¹³C-octanoic acid breath test: a noninvasive technique to assess gastric emptying in preterm infants. *J Pediatr Gastroenterol Nutr.* 1996 Aug;23(2):111-7.
19. Ziegler D, Schadewaldt P, Pour Mirza A, Piolot R, Schommartz B, Reinhardt M, Vosberg H, Brosicke H, Gries FA. [¹³C]octanoic acid breath test for non-invasive assessment of gastric emptying in diabetic patients: validation and relationship to gastric symptoms and cardiovascular autonomic function. *Diabetologia.* 1996 Jul;39(7):823-30.
20. Maes BD, Ghoo YF, Rutgeerts PJ, Hiele MI, Geypens B, Vantrappen G. [¹³C]octanoic acid breath test to measure gastric emptying rate of solids. *Dig Dis Sci.* 1994 Dec;39(12 Suppl):104S-106S.
21. Maes BD, Ghoo YF, Geypens BJ, Mys G, Hiele MI, Rutgeerts PJ, Vantrappen G. Combined carbon-¹³-glycine/carbon-14-octanoic acid breath test To monitor gastric emptying rates of liquids and solids. *J Nucl Med.* 1994 May;35(5):824-31.
22. Maes BD, Hiele MI, Geypens BJ, Rutgeerts PJ, Ghoo YF, Vantrappen G. Pharmacological modulation of gastric emptying rate of solids as measured by the carbon labelled octanoic acid breath test: influence of erythromycin and propantheline. *Gut.* 1994 Mar;35(3):333-7.
23. Ghoo YF, Maes BD, Geypens BJ, Mys G, Hiele MI, Rutgeerts PJ, Vantrappen G. Measurement of gastric emptying rate of solids by means of a carbon-labeled octanoic acid breath test. *Gastroenterol.* 1993 Jun;104(6):1640-7.

3. Gastric Emptying of Liquids ¹³C-Acetate Breath Test

■ Recommended literature

1. Shimamoto C, Hirata I, Hiraike Y, Takeuchi N, Nomura T, Katsu K. Evaluation of gastric motor activity in the elderly by electrogastrography and the ¹³C-Acetate breath test. *Gerontology*. 2002 Nov-Dec;48(6):381-6.
2. Urita Y, Hike K, Torii N, Kikuchi Y, Sasajima M, Miki K. Efficacy of lactulose plus ¹³C-acetate breath test in the diagnosis of gastrointestinal motility disorders. *J Gastroenterol*. 2002;37(6):442-8
3. Gonzalez A, Mugueta C, Parra D, Labayen I, Martinez A, Varo N, Monreal I, Gil MJ. Characterisation with stable isotopes of the presence of a lag phase in the gastric emptying of liquids. *Eur J Nutr*. 2000 Oct;39(5):224-8.
4. van Nieuwenhoven MA, Wagenmakers AJ, Senden JM, Brouns F, Brummer RJ. Performance of the [¹³C]-acetate gastric emptying breath test during physical exercise. *Eur J Clin Invest*. 1999 Nov;29(11):922-8.
5. Mudambo KS, Leese GP, Rennie MJ. Gastric emptying in soldiers during and after field exercise in the heat measured with the [¹³C]acetate breath test method. *Eur J Appl Physiol Occup Physiol*. 1997;75(2):109-14.
6. Braden B, Adams S, Duan LP, Orth KH, Maul FD, Lembcke B, Hor G, Caspary WF. The [¹³C]acetate breath test accurately reflects gastric emptying of liquids in both liquid and semisolid test meals. *Gastroenterology*. 1995 Apr;108(4):1048-55.
7. Leese GP, Bowtell J, Mudambo S, Reynolds N, Thompson J, Srimgeour CM, Rennie MJ. Post-exercise gastric emptying of carbohydrate solutions determined using the ¹³C-Acetate breath test. *Eur J Appl Physiol Occup Physiol*. 1995;71(4):306-10.
8. Mossi S, Meyer-Wyss B, Beglinger C, Schwizer W, Fried M, Ajami A, Brignoli R. Gastric emptying of liquid meals measured noninvasively in humans with ¹³C-Acetate breath test. *Dig Dis Sci*. 1994 Dec;39(12 Suppl):107S-109S.