



# Analysis of Per- and Polyfluoroalkyl Substances in Drinking Water Using Sequential and Parallel Automated Solid Phase Extraction Using EPA Method 537.1

Ruud Addink and Nem Kaurin  
Fluid Management Systems  
Watertown MA



- Per- and polyfluoralkyl substances (PFAS) contain a perfluorinated or polyfluorinated carbon chain moiety such as  $F(CF_2)_n-$  or  $F(CF_2)_n-(C_2H_4)_n$ .
- These make up a large group of persistent chemicals used in industrial processes and consumer applications:
  - Stain-Resistant Coatings for textiles and carpets
  - Grease-Proof Coatings for paper products approved for food contact
  - Firefighting Foams
  - Mining and Oil Well Surfactants
  - Floor Polishes
  - Insecticide formulations

# Origin


- Industrial Sites
  - Airport Fire Training Areas
  - Wastewater Treatment Facilities
  - Widespread use for over 60 years
  - Very resistant to degradation
  - Ubiquitous Compound in the Environment
- 

- **Human exposure is linked to adverse effects**
  - Developmental issues in off-spring
  - Cancer
  - Immune system suppression
  - Endocrine disruption
  - Elevated levels of Cholesterol
  - Obesity


# Source concerns

- Many water sources worldwide are found to be contaminated.
- Two compounds most studied:
  - Perfluorooctane sulphonate (PFOS)
  - Perfluorooctanoic acid (PFOA)
- Millions have been exposed through drinking water supplies. EPA advisory levels < 70ng/L

# Regulation

- PFAS regulation is emerging worldwide
  - Stockholm Convention review has recommended (partial) elimination
  - EPA has not issued Maximum Contaminant Levels (MCLs)
- 

# The Analysis of PFAS

- Many of Thousands Samples are now being analyzed and more locations are starting to be analyzed for PFAS
    - Drinking Water
    - Waste Water
    - Human Serum
    - Biota
    - Soils
- 

# Challenges in the Analysis of PFAS

- The Analytical Systems are expensive
  - UPLC/MS systems
    - Require expertise in a new technology
- Manual Sample Prep processes
  - Inconsistent results
  - Elevated Background issues
  - Labor intensive
  - Extraction can take up to 2 hours
  - Concentration can take up to 2 hours



- Automate the Sample Prep Workflow
  - Automate the Solid Phase Extraction Step
  - Automate the Concentration/Evaporation Step
- Automated SPE extractions and Concentration is a very green technique
  - Reduces Solvent Use
  - Reduces Solvent Disposal Costs
  - Reduces Solvent emissions
- FMS automated SPE systems deliver consistent, reproducible results
- Solid Phase Extraction is a well accepted technology

# Reasons for SPE

- Reduced solvent
- Reduced glassware
- Simplified faster procedures (80 min automated vs 150 min manual)
- Automation versus manual protocols = Reproducibility



# Determining Factors

- Ability to load samples by both positive pressure and vacuum.
- Ability to dry cartridges by both vacuum and positive gas pressure (N<sub>2</sub>).
- Easily handle a wide variety of cartridge designs and sizes without cumbersome modifications.

# Automated SPE System for PFAS extraction (1)



- Expandable from 1 to 6 modules
- Parallel and Sequential Extraction
- Direct to Concentrator and Vial
- All Inert Peek and Stainless Steel Surfaces

# Automated SPE System for PFAS extraction (2)

- Low Background system
  - Peek and Stainless components
- Modular and Expandable System
  - Up to 6 modules
- High Throughput Runs Sample Extraction in Parallel and Sequential mode
- Up to 30 samples run unattended in 6 h period


# Automated SPE System for PFAS extraction (3)

- Uses Vacuum for Sample Loading
- Uses Positive Pressure Pumping for Precise delivery of Elution and Wash Solvent





**No Teflon**




# TurboTrace PFAS







# EZPFC™ (1)



# EZPFC™ (2)



# Extraction procedure (1)

- 250 mL water samples were spiked with 2, 50 or 100 ng/L PFAS standards
  - Use FMS 500 mg cartridge (SDVB).
  - Condition cartridge with 15 mL methanol (soak 2 min).
  - Condition cartridge with 18 mL water (soak 2 min).
- 

## Extraction procedure (2)

- Load samples onto system
- Load samples across cartridges under -8 inches Hg vacuum (20-25 min)
- Rinse bottle with 2 x 7.5 mL of water and load onto the cartridge under negative pressure.

## Extraction procedure (3)

- Dry cartridges under nitrogen until no residual water is present (5 min)
- Rinse the sample bottles and elute with 2 x 4 mL methanol




- Pre-heat temp: 60-65 °C
- Pre-heat time: 20 minutes
- Heat in Time mode: 60-65 °C
- Nitrogen pressure: 8 PSI
- The extracts were concentrated to dryness, reconstituted with methanol:water, add internal standard. The samples had a final volume of 1 mL for LC/MS analysis.

- **UPLC Conditions**


- Waters Acquity UPLC with Q-TOF (Xevo G2-XS) and HR-MS
- Acquity HSS T3 column (2.1 mm × 100 mm, 1.8 μm)
- Negative ESI

# Recoveries 2 ppt PFAS







# Recoveries 50 ppt PFAS



# Recoveries 100 ppt PFAS



# PFAS Background



# Conclusions

- It is possible to automate the sample preparation of Per- and Polyfluoroalkyl Substances with the FMS SPE Parallel/Sequential System and SuperVap Concentrator for high throughput analysis
- Delivers consistent and reproducible results for PFAS analysis
- Can run up to 30 samples fully automated and unattended over a 6 h period
- The system, by design, has very low background PFAS allowing for analysis of samples without any significant interference.
- All models of FMS SPE systems are available as PFAS systems

# Questions

- Questions?
- See us at booth # 4043.

