



NEW PRODUCTS

U.S. EPA Method 537 Analyte Primary Dilution Standard Solution/Mixtures & Linear/Branched Solution/Mixtures of N-MeFOSAA and N-EtFOSAA

The U.S. Environmental Protection Agency developed Method 537 (Version 1.1) for the determination of selected perfluorinated alkyl acids in drinking water by solid phase extraction and liquid chromatography/tandem mass spectrometry (LC/MS/MS). This method requires the use of branched/linear isomer mixtures of PFHxS, PFOS, N-MeFOSAA, and N-EtFOSAA. However, until now, characterized branched/linear mixtures of N-MeFOSAA and N-EtFOSAA were not commercially available. Therefore, in response to market demand, **Wellington** is pleased to announce the development and release of linear/branched isomer mixtures of N-MeFOSAA (**br-NMeFOSAA**) and N-EtFOSAA (**br-NEtFOSAA**) which have been characterized as to their isomeric content by ¹⁹F NMR.

To facilitate the use of these new branched/linear standards according to EPA requirements, **Wellington** has also prepared two analyte primary dilution standard solutions for U.S. EPA Method 537. **EPA-537PDS** contains all of the native PFAS analytes required by Method 537 with PFHxS, PFOS, N-MeFOSAA, and N-EtFOSAA being present as linear/branched mixtures whereas **EPA-537PDS-L** contains only linear isomers of the same components.

	Catalogue Number	Product (methanol)	Qty	Conc
NEW	br-NMeFOSAA	N-Methylperfluorooctanesulfonamidoacetic acid Isomeric Mixture	1.2 ml	50 µg/ml
NEW	br-NEtFOSAA	N-Ethylperfluorooctanesulfonamidoacetic acid Isomeric Mixture	1.2 ml	50 µg/ml

	Catalogue Number	Product (methanol)	Qty	Conc
NEW	EPA-537PDS	Analyte Primary Dilution Standard (branched/linear mix)	1.2 ml	2.0 µg/ml ea
NEW	EPA-537PDS-L	Analyte Primary Dilution Standard (linear isomers only)	1.2 ml	2.0 µg/ml ea
	EPA-537SS	Surrogate Primary Dilution Standard (SUR PDS)	1.2 ml	
	MPFHxA	Perfluoro-n-[1,2- ¹³ C ₂]hexanoic acid		1.0 µg/ml
	MPFDA	Perfluoro-n-[1,2- ¹³ C ₂]decanoic acid		1.0 µg/ml
	d5-N-EtFOSAA	N-ethyl-d ₅ -perfluoro-1-octanesulfonamidoacetic acid		4.0 µg/ml
	EPA-537IS	Internal Standard Primary Dilution Standard (IS PDS)	1.2 ml	
	M2PFOA	Perfluoro-n-[1,2- ¹³ C ₂]octanoic acid		1.0 µg/ml
	MPFOS	Sodium perfluoro-1-[1,2,3,4- ¹³ C ₄]octanesulfonate		3.0 µg/ml
	d3-N-MeFOSAA	N-methyl-d ₃ -perfluoro-1-octanesulfonamidoacetic acid		4.0 µg/ml

Please contact your local distributor or info@well-labs.com for pricing and delivery.

Visit our website (www.well-labs.com) for a complete listing of our new products.



LC/MS/MS Chromatograms

